SINEKRON-S

PANESTETIS

POWe

Infting

CORDLESS ULTRASONIC SPATULA



Ideal for a perfect peeling, deep cleansing and delivery of product.



SUPERFICIAL EXFOLIATION **AND PEELING**

> For a pre-treatment face cleansing.



DEEP FACE CLEANSING

Thanks to vibrations which do not cause stress or damage the tissue.



DEEP DELIVERY **OF THE PRODUCT**

Delivery of active ingredients for an even more effective result.

Perfect with

LIFTING MASK

Pure biocellulose fiber mask which provides an immediate lifting effect.

It adheres perfectly to the face and releases active ingredients with lifting and anti-aging properties to restore tone, firmness hydration and brightness already from the first application.

Ideal for bidirectional action as a final step after the treatment carried out with technologies.

ELE@TRA STANS MORE ADDRESS AND

BRILLIANT MASK

Pure biocellulose fiber mask yhat ensures a professional treatment to correct hyperpigmentation problems.

Its formula includes lightening active ingredients designed to reduce and minimize localized spots in the epidermis.

It brightens, illuminates and evens complexion.



SINEKRON-S **CORDLESS ULTRASONIC SPATULA**



SONIC: ACTIVATES THE CONTINUOUS

- LIFTING: ACTIVATES THE CONTINUOUS ULTRASOUND

- ION +: ACTIVATES THE PULSED ULTRASOUND MODE AND THE POSITIVE (+) OR NEGATIVE (-) IONIZATION

LOW HIGH: ENABLES TO ADJUST THE INTENSITY OF THE ULTRASOUND (LOW OR HIGH) THIS FUNCTION IS AVAILABLE FOR ALL THE PROGRAMS



Panestetic, certified ISO 9001 and ISO 13485, is a leading Italian manufacturer of aesthetic and medical devices since 1984.

Thanks to the great experience gained in the research and development of new technologies and devices, the company has developed a range of professional cosmetic products to work in synergy with beauty machines, boosting the effects in terms of performances and more visible results.

PANESTETIC SRL

Via Vecchia Francesca 10/A – 24040 Ciserano (BG) Italy ph. +39 035 883087 | www.panestetic.com | export@panestetic.it PANESTETICTECHNOLOGY @PANESTETIC